

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

EPT 1 Conductivity cell - Epoxy body - for difficult samples

Item Code: 50003382



DATA SHEET

Features

Cable	Si
Connection	BNC / RCA
	Epoxy C=1 cell, carbon rings, with built-in temperature sensor.
Description	Fixed cable 1 m, BNC-Cinch. Range 0...100 mS
Diaphragm	Carbon
Body diameter	12 mm
Image 2	https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6
Body lenght	120 mm
Cable length	1 m
Body material	Epoxy
Name	EPT 1
Cond measuring range	10 µS/cm ... 100 mS/cm
temperature measurement range	0 ... +60 °C
Temperature sensor	Si / spina RCA Cinch