

Electrode RX 3mt cable for REVIO and PCD 650

Item Code: 51004132



DATA SHEET

Features

Features	Maintenance free
Connection	BNC
Description	Combined Redox electrode. General purpose plastic body. 3 m fixed cable with BNC plug. Gel filling with porous septum. Maintenance-free. Working range: $\pm 1,000$ mV, temperature 0...60 °C
Diaphragm	porous septum
Body diameter	12 mm
Electrolyte	Gel
Image 2	Redox electrode
Body lenght	120 mm
Cable length	3 m
Body material	Plastic
Name	Redox electrode
mV ORP measuring range	-1000 ... +1000 mV